

DICOM Conformance Statement

RETI-Port/Scan 21



ROLAND CONSULT
Electrophysiology and Imaging

No. 01
Rev. 01

DICOM Conformance Statement

RETI-Port/Scan 21

Number 01

Revision 01



	Name	Bereich	Datum	Unterschrift
Erstellt:	S. Schmidt	E	15.03.2021	
Geprüft/ Genehmigt:	J. Finger	QB	15.03.2021	

1 CONFORMANCE STATEMENT OVERVIEW

The RETI-Port/Scan 21 implements the necessary DICOM services to download work lists from an information system, convert PDF documents with measured values, curves and derivatives into the DICOM format "Encapsulated PDF" and transfer them to a DICOM-enabled network storage device (typically a PACS).

Table B.1-1 provides an overview of the network services supported by RETI-Port/Scan 21.

Table B.1-1
NETWORK SERVICES

SOP Classes	Use of Service (SCU)	Provider of Service (SCP)
Verification	Yes	Yes
Transfer		
Encapsulated PDF Storage	Yes	No
Workflow Management		
Modality Worklist	Yes	No

DICOM Conformance Statement

RETI-Port/Scan 21



No. 01
Rev. 01

2 Table of Contents

1	CONFORMANCE STATEMENT OVERVIEW	2
2	Table of Contents	3
	List of Figures	4
	List of Tables	4
3	INTRODUCTION	6
3.1	AUDIENCE	6
3.2	REMARKS	6
3.3	TERMS AND DEFINITIONS	6
3.4	BASICS OF DICOM COMMUNICATION	8
3.5	ABBREVIATIONS, REFERENCES	8
3.6	REFERENCES	10
4	NETWORKING	10
4.1	IMPLEMENTATION MODEL	10
4.1.1	Application Data Flow	10
4.2	Functional Definition of AE's	11
4.2.1	Functional Definition of Verification Application Entity	11
4.2.2	Functional Definition of Send Application Entity	11
4.2.3	Functional Definition of Workflow Application Entity	11
4.2.4	Sequencing of Real-World Activities	12
4.3	AE SPECIFICATIONS	12
4.3.1	Send Application Entity Specification	12
4.3.1.1	SOP Classes	12
4.3.1.2	Association Policies	13
4.3.1.2.1	General	13
4.3.1.2.2	Number of Associations	13
4.3.1.2.3	Asynchronous Nature	13
4.3.1.2.4	Implementation Identifying Information	13
4.3.1.3	Association Initiation Policy	13
4.3.1.3.1	Activity – Send Images	13
4.3.2	Workflow Application Entity Specification	15
4.3.2.1	SOP Classes	15
4.3.2.2	Association Policies	16
4.3.2.2.1	General	16
4.3.2.2.2	Number of Associations	16
4.3.2.2.3	Asynchronous Nature	16
4.3.2.2.4	Implementation Identifying Information	16
4.3.2.1	Association Initiation Policy	17
4.3.2.1.1	Activity – Worklist Update	17
4.3.3	Verification Application Entity	23
4.3.3.1	SOP Classes	23
4.3.3.2	Association Policies	23
4.3.3.2.1	General	23
4.3.3.2.2	Number of Associations	23
4.3.3.2.3	Asynchronous Nature	23
4.3.3.2.4	Implementation Identifying Information	23
4.3.3.3	Association Initiation and Acceptance Policy	24
4.3.3.3.1	Activity: Verify connectivity with remote AE	24
4.4	NETWORK INTERFACES	25
4.4.1	Physical Network Interface	25
4.4.2	Additional Protocols	25
5	CONFIGURATION	26
5.1	AE Title/Presentation Address Mapping	26

No. 01
Rev. 01

5.1.1	Local AE Titles	26
5.2	Parameters.....	26
6	SUPPORT OF CHARACTER SETS	28
7	SECURITY	29
8	ANNEXES.....	30
8.1	IOD CONTENTS	30
8.1.1	Created SOP Instances	30
8.1.1.1	Encapsulated PDF IOD.....	30
8.1.2	IOD Module Definitions	31
8.1.2.1	Common Modules	31
8.1.2.1	Secondary Capture Image specific Modules	33
8.1.2.2	Encapsulated PDF Specific Attributes	33
8.1.4	Attribute mapping	34
8.1.5	Coerced/Modified Fields	35
8.2	PRIVATE TRANSFER SYNTAXES	36
8.3	Association Rejection Reasons	37
9	Revision History	37

List of Figures

Figure 4-1	RETI-Port/Scan 21 DATA FLOW DIAGRAM.....	10
Figure 4-2	SEQUENCE DIAGRAM – IMAGE ACQUISITION AND ARCHIVING TRIGGERED BY WORKLIST 12	
Figure 4-3	SEQUENCING OF ACTIVITY – WORKLIST UPDATE	17
Figure 4-4	SEQUENCING OF ACTIVITY: VERIFY CONNECTIVITY WITH REMOTE AE	24

List of Tables

Table 4.3-1	SOP CLASS FOR SEND AE	12
Table 4.3-2	DICOM APPLICATION CONTEXT FOR SEND AE.....	13
Table 4.3-3	NUMBER OF ASSOCIATIONS INITIATED FOR SEND AE	13
Table 4.3-4	ASYNCHRONOUS NATURE AS A SCU FOR SEND AE	13
Table 4.3-5	DICOM IMPLEMENTATION CLASS AND VERSION FOR SEND AE	13
Table 4.3-6	PROPOSED PRESENTATION CONTEXTS FOR ACTIVITY SEND IMAGES	14
Table 4.3-7	STORAGE C-STORE RESPONSE STATUS HANDLING BEHAVIOR.....	14
Table 4.3-8	STORAGE COMMUNICATION FAILURE BEHAVIOR	15
Table 4.3-9	SOP CLASSES FOR WORKFLOW AE	16
Table 4.3-10	DICOM APPLICATION CONTEXT FOR WORKFLOW AE	16
Table 4.3-11	NUMBER OF ASSOCIATIONS INITIATED FOR WORKFLOW AE	16
Table 4.3-12	ASYNCHRONOUS NATURE AS A SCU FOR WORKFLOW AE	16
Table 4.3-13	DICOM IMPLEMENTATION CLASS AND VERSION FOR WORKFLOW AE	16
Table 4.3-14	PROPOSED PRESENTATION CONTEXTS FOR ACTIVITY WORKLIST UPDATE	18
Table 4.3-15	MODALITY WORKLIST - FIND SOP Class - C-FIND REQUEST IDENTIFIER	19
Table 4.3-16	MODALITY WORKLIST C-FIND RESPONSE STATUS HANDLING BEHAVIOR.....	22
Table 4.3-17	MODALITY WORKLIST COMMUNICATION FAILURE BEHAVIOR.....	22
Table 4.3-18	DICOM APPLICATION CONTEXT FOR VERIFICATION AE	23
Table 4.3-19	NUMBER OF SIMULTANEOUS ASSOCIATIONS FOR VERIFICATION AE	23
Table 4.3-20	ASYNCHRONOUS NATURE AS A SCU FOR VERIFICATION AE	23
Table 4.3-21	DICOM IMPLEMENTATION CLASS AND VERSION FOR VERIFICATION AE	24
Table 4.3-22	PROPOSED PRESENTATION CONTEXTS FOR VERIFICATION AE	25
Table 5.1-1	AE TITLE CONFIGURATION TABLE	26
Table 5.2-1	CONFIGURATION PARAMETERS TABLE	26
Table 8.1-1	IOD OF CREATED ENCAPSULATED PDF SOP INSTANCES	30
Table 8.1-2	PATIENT MODULE ATTRIBUTES OF CREATED SOP INSTANCES	31
Table 8.1-3	GENERAL STUDY MODULE ATTRIBUTES OF CREATED SOP INSTANCES	31
Table 8.1-4	PATIENT STUDY MODULE ATTRIBUTES OF CREATED SOP INSTANCES	32

DICOM Conformance Statement

RETI-Port/Scan 21



ROLAND CONSULT
Electrophysiology and Imaging

No. 01

Rev. 01

Table 8.1-5 GENERAL SERIES MODULE ATTRIBUTES OF CREATED SOP INSTANCES	32
Table 8.1-6 GENERAL EQUIPMENT MODULE ATTRIBUTES OF CREATED SOP INSTANCES	32
Table 8.1-7 SOP COMMON MODULE ATTRIBUTES OF CREATED SOP INSTANCES	33
Table 8.1-8 SC EQUIPMENT MODULE ATTRIBUTES OF CREATED SOP INSTANCES	33
Table 8.1-9 ENCAPSULATED DOCUMENT SERIES MODULE ATTRIBUTES OF CREATED SOP INSTANCES	33
Table 8.1-10 ENAPSULATED DOCUMENT MODULE ATTRIBUTES OF CREATED SOP INSTANCES	34
Table 8.1-11 ATTRIBUTE MAPPING BETWEEN MODALITY WORKLIST AND CREATED INSTANCES	34
Table 8.3-1 ASSOCIATION REJECTION REASONS	37

 ROLAND CONSULT Electrophysiology and Imaging	DICOM Conformance Statement RETI-Port/Scan 21
No. 01 Rev. 01	

3 INTRODUCTION

3.1 AUDIENCE

This document is intended for people that need to understand how RETI-Port/Scan 21 will integrate into their healthcare facility. This includes both those responsible for overall imaging network policy and architecture, as well as integrators who need to have a detailed understanding of the DICOM features of the product. This document contains some basic DICOM definitions so that any reader may understand how this product implements DICOM features. However, integrators are expected to fully understand all the DICOM terminology, how the tables in this document relate to the product's functionality, and how that functionality integrates with other devices that support compatible DICOM features.

3.2 REMARKS

The scope of this DICOM Conformance Statement is to facilitate integration between RETI-Port/Scan 21 and other DICOM products. The Conformance Statement should be read and understood in conjunction with the DICOM Standard. DICOM by itself does not guarantee interoperability. However, the Conformance Statement facilitates a first-level comparison for interoperability between different applications supporting compatible DICOM functionality.

This Conformance Statement is not supposed to replace validation with other DICOM equipment to ensure proper exchange of intended information. In fact, the user should be aware of the following important issues:

- The comparison of different Conformance Statements is just the first step towards assessing interconnectivity and interoperability between the product and other DICOM conformant equipment.
- Test procedures should be defined and executed to validate the required level of interoperability with specific compatible DICOM equipment, as established by the healthcare facility.

3.3 TERMS AND DEFINITIONS

Terms and definitions should be listed here. The following example may be used as a template:

Informal definitions are provided for the following terms used in this Conformance Statement. The DICOM Standard is the authoritative source for formal definitions of these terms.

Abstract Syntax	The information agreed to be exchanged between applications, generally equivalent to a Service/Object Pair (SOP) Class. Examples: Verification SOP Class, Modality Worklist Information Model Find SOP Class, Computed Radiography Image Storage SOP Class.
Application Entity (AE)	An end point of a DICOM information exchange, including the DICOM network or media interface software, i.e., the software that sends or receives DICOM information objects or messages. A single device may have multiple Application Entities.
Application Entity Title (AET)	The externally known name of an Application Entity, used to identify a DICOM application to other DICOM applications on the network.
Application Context	The specification of the type of communication used between Application Entities. Example: DICOM network protocol.
Association	A network communication channel set up between Application Entities.
Attribute	A unit of information in an object definition; a data element identified by a tag. The information may be a complex data structure (Sequence), itself composed of lower-level data elements. Examples: Patient ID (0010,0020), Accession Number (0008,0050), Photometric Interpretation (0028,0004), Procedure Code Sequence (0008,1032).

DICOM Conformance Statement

RETI-Port/Scan 21



No. 01
Rev. 01

Information Object Definition (IOD)	The specified set of Attributes that comprise a type of data object; does not represent a specific instance of the data object, but rather a class of similar data objects that have the same properties. The Attributes may be specified as Mandatory (Type 1), Required but possibly unknown (Type 2), or Optional (Type 3), and there may be conditions associated with the use of an Attribute (Types 1C and 2C). Examples: MR Image IOD, CT Image IOD, Print Job IOD.
Joint Photographic Experts Group (JPEG)	A set of standardized image compression techniques, available for use by DICOM applications.
Media Application Profile	The specification of DICOM information objects and encoding exchanged on removable media (e.g., CDs).
Module	A set of Attributes within an Information Object Definition that are logically related to each other. Example: Patient Module includes Patient Name, Patient ID, Patient Birth Date, and Patient Sex.
Negotiation	First phase of Association establishment that allows Application Entities to agree on the types of data to be exchanged and how that data will be encoded.
Presentation Context	The set of DICOM network services used over an Association, as negotiated between Application Entities; includes Abstract Syntaxes and Transfer Syntaxes.
Protocol Data Unit (PDU)	A packet (piece) of a DICOM message sent across the network. Devices must specify the maximum size packet they can receive for DICOM messages.
Security Profile	A set of mechanisms, such as encryption, user authentication, or digital signatures, used by an Application Entity to ensure confidentiality, integrity, and/or availability of exchanged DICOM data
Service Class Provider (SCP)	Role of an Application Entity that provides a DICOM network service; typically, a server that performs operations requested by another Application Entity (Service Class User). Examples: Picture Archiving and Communication System (image storage SCP, and image query/retrieve SCP), Radiology Information System (modality worklist SCP).
Service Class User (SCU)	Role of an Application Entity that uses a DICOM network service; typically, a client. Examples: imaging modality (image storage SCU, and modality worklist SCU), imaging workstation (image query/retrieve SCU)
Service/Object Pair Class (SOP Class)	The specification of the network or media transfer (service) of a particular type of data (object); the fundamental unit of DICOM interoperability specification. Examples: Ultrasound Image Storage Service, Basic Grayscale Print Management.
Service/Object Pair Instance (SOP Instance)	An information object; a specific occurrence of information exchanged in a SOP Class. Examples: a specific x-ray image.
Tag	A 32-bit identifier for a data element, represented as a pair of four-digit hexadecimal numbers, the "group" and the "element". If the "group" number is odd, the tag is for a private (manufacturer-specific) data element. Examples: (0010,0020) [Patient ID], (07FE,0010) [Pixel Data], (0019,0210) [private data element]

 ROLAND CONSULT Electrophysiology and Imaging	DICOM Conformance Statement RETI-Port/Scan 21
No. 01 Rev. 01	

Transfer Syntax	The encoding used for exchange of DICOM information objects and messages. Examples: JPEG compressed (images), little endian explicit value representation.
Unique Identifier (UID)	A globally unique "dotted decimal" string that identifies a specific object or a class of objects; an ISO-8824 Object Identifier. Examples: Study Instance UID, SOP Class UID, SOP Instance UID.
Value Representation (VR)	The format type of an individual DICOM data element, such as text, an integer, a person's name, or a code. DICOM information objects can be transmitted with either explicit identification of the type of each data element (Explicit VR), or without explicit identification (Implicit VR); with Implicit VR, the receiving application must use a DICOM data dictionary to look up the format of each data element.

3.4 BASICS OF DICOM COMMUNICATION

This section describes terminology used in this Conformance Statement for the non-specialist. This section is not a substitute for training about DICOM, and it makes many simplifications about the meanings of DICOM terms.

Two Application Entities (devices) that want to communicate with each other over a network using DICOM protocol must first agree on several things during an initial network "handshake". One of the two devices must initiate an Association (a connection to the other device), and ask if specific services, information, and encoding can be supported by the other device (Negotiation).

DICOM specifies several network services and types of information objects, each of which is called an Abstract Syntax for the Negotiation. DICOM also specifies a variety of methods for encoding data, denoted Transfer Syntaxes. The Negotiation allows the initiating Application Entity to propose combinations of Abstract Syntax and Transfer Syntax to be used on the Association; these combinations are called Presentation Contexts. The receiving Application Entity accepts the Presentation Contexts it supports.

For each Presentation Context, the Association Negotiation also allows the devices to agree on Roles - which one is the Service Class User (SCU - client) and which is the Service Class Provider (SCP - server). Normally the device initiating the connection is the SCU, i.e., the client system calls the server, but not always.

The Association Negotiation finally enables exchange of maximum network packet (PDU) size, security information, and network service options (called Extended Negotiation information).

The Application Entities, having negotiated the Association parameters, may now commence exchanging data. Common data exchanges include queries for worklists and lists of stored images, transfer of image objects and analyses (structured reports), and sending images to film printers. Each exchangeable unit of data is formatted by the sender in accordance with the appropriate Information Object Definition, and sent using the negotiated Transfer Syntax. There is a Default Transfer Syntax that all systems must accept, but it may not be the most efficient for some use cases. Each transfer is explicitly acknowledged by the receiver with a Response Status indicating success, failure, or that query or retrieve operations are still in process.

Two Application Entities may also communicate with each other by exchanging media (such as a CD-R). Since there is no Association Negotiation possible, they both use a Media Application Profile that specifies "pre-negotiated" exchange media format, Abstract Syntax, and Transfer Syntax.

3.5 ABBREVIATIONS, REFERENCES

Definitions, terms, and abbreviations used in this document are defined within the different parts of the DICOM standard.

Additional Abbreviations and terms are as follows:

DICOM Conformance Statement

RETI-Port/Scan 21



No. 01
Rev. 01

AE	DICOM Application Entity
AET	Application Entity Title
ASCII	American Standard Code for Information Interchange
CSE	Customer Service Engineer
DB	Database
DCS	DICOM Conformance Statement
DICOM	Digital Imaging and Communications in Medicine
FSC	File Set Creator
FSR	File Set Reader
FSU	File Set Updater
GSDF	Grayscale Standard Display Function
HIS	Hospital Information System
IHE	Integrating the Healthcare Enterprise
IOD	DICOM Information Object Definition
ISO	International Standard Organization
MPPS	Modality Performed Procedure Step
n. a.	not applicable
NEMA	National Electrical Manufacturers Association
PACS	Picture Archiving and Communication System
PDU	DICOM Protocol Data Unit
RIS	Radiology Information System
RP	Reference Point
SC	Secondary Capture
SCU	DICOM Service Class User (DICOM client)
SCP	DICOM Service Class Provider (DICOM Server)
SOP	DICOM Service-Object Pair
SPS	Scheduled Procedure Step
SR	Structured Report
TFT	Thin Film Transistor (Display)
TID	Template ID
UID	Unique Identifier
UTF-8	Unicode Transformation Format-8
VR	Value Representation

 ROLAND CONSULT Electrophysiology and Imaging	DICOM Conformance Statement RETI-Port/Scan 21
No. 01 Rev. 01	

3.6 REFERENCES

[NEMA PS3] Digital Imaging and Communications in Medicine (DICOM) Standard, available free at <http://medical.nema.org/>¹

4 NETWORKING

4.1 IMPLEMENTATION MODEL

4.1.1 Application Data Flow




Figure 4-1 RETI-Port/Scan 21 DATA FLOW DIAGRAM

Conceptually the network services may be modeled as the following separate AEs, though in fact all the AEs share a single (configurable) AE Title:

The Send Application Entity sends Encapsulated PDF objects to a remote AE. It is associated with the real-world activity "Send Images" which is performed after completion of an exam.

The Workflow Application Entity receives worklist information from a remote AE. It is associated with the real-world activity "Update Worklist". Update Worklist is performed as a result of an operator request. The Workflow AE queries a remote AE for worklist items matching a query request defined by the operator (time range of scheduled examination, PACS id of patient or * for all open WLEs).

The Verification AE sends verification requests to and answers verification requests from a remote AE.

¹ The DICOM Standard is under continuous maintenance, the current official version is available at <http://dicom.nema.org>

DICOM Conformance Statement

RETI-Port/Scan 21



No. 01
Rev. 01

4.2 Functional Definition of AE's

4.2.1 Functional Definition of Verification Application Entity

The DICOM Verification AE is used to check if the configuration of the remote AE is correct, and the remote AE is listening. The user can trigger the echo request. The Request is initiated in the DICOM window ("verify" button)

The result of the request is displayed immediately. Furthermore, the DICOM Verification AE responds to verification requests from remote AEs. Technically, this part is included in the system service, which is also hosting the Storage Service AE.

Technically, the Verification SCU is an integral component of all SCUs inside the RETI-Port/Scan 21, so every SCU can verify basic connectivity before putting a Service Class specific request towards the SCP.

4.2.2 Functional Definition of Send Application Entity

If a send-job is triggered for a file or a study, an association request is sent to the destination AE and upon successful negotiation of a Presentation Context, the document transfer is started. If the association cannot be opened or the transfer of the acquired documents fails for any other reason, the failed transfer attempt will be recorded in RETI-Port/Scan 21's database. A configurable number of retries will be executed in configurable time intervals until either the transfer succeeds or the maximum number of retries is reached.

The user can display the transfer protocol.

In case of failed transfers, a notification is displayed when the DICOM window is launched.

4.2.3 Functional Definition of Workflow Application Entity

The Workflow Application Entity attempts to download a Worklist from a remote node. If the Workflow AE establishes an Association to a remote AE, it will transfer all worklist items via the open Association. The results will be displayed in a separate list, which will be cleared with the next Worklist Update.

4.2.4 Sequencing of Real-World Activities




Figure 4-2 SEQUENCE DIAGRAM – IMAGE ACQUISITION AND ARCHIVING TRIGGERED BY WORKLIST

Figure 4-2 shows the sequence of activities to receive worklist items, create images and archive them. Other workflow situations (e.g., unscheduled procedure steps) will have other sequencing constraints.

4.3 AE SPECIFICATIONS

4.3.1 Send Application Entity Specification

4.3.1.1 SOP Classes

The Send Application Entity provides Standard Conformance to the following SOP Class:

Table 4.3-1 SOP CLASS FOR SEND AE

SOP Class Name	SOP Class UID	SCU	SCP
Encapsulated PDF Storage	1.2.840.10008.5.1.4.1.1.104.1	Yes	No

DICOM Conformance Statement

RETI-Port/Scan 21



No. 01
Rev. 01

4.3.1.2 Association Policies

4.3.1.2.1 General

The DICOM standard application context name for DICOM 3.0 is always proposed:

Table 4.3-2 DICOM APPLICATION CONTEXT FOR SEND AE

Application Context Name	1.2.840.10008.3.1.1.1
--------------------------	-----------------------

4.3.1.2.2 Number of Associations

The Send AE initiates one Association at a time for each destination to which a transfer request is being processed in the active job queue list. Only one job will be active at a time, the other remains pending until the active job is completed or failed.

Table 4.3-3 NUMBER OF ASSOCIATIONS INITIATED FOR SEND AE

Maximum number of simultaneous Associations	1
---	---

4.3.1.2.3 Asynchronous Nature

The Send AE does not support asynchronous communication (multiple outstanding transactions over a single Association).

Table 4.3-4 ASYNCHRONOUS NATURE AS A SCU FOR SEND AE

Maximum number of outstanding asynchronous transactions	1
---	---

4.3.1.2.4 Implementation Identifying Information

The implementation information for this Application Entity is:

Table 4.3-5 DICOM IMPLEMENTATION CLASS AND VERSION FOR SEND AE

Implementation Class UID	1.2.826.0.1.3680043.2.891.113
Implementation Version Name	DICOM_CONNECT_25

4.3.1.3 Association Initiation Policy

4.3.1.3.1 Activity – Send Images

4.3.1.3.1.1 Description and Sequencing of Activities

A user can select image and non-image objects and request them to be sent to a configured destination. SOP instances of a specific study can be sent together or in different associations.

If a send-job is triggered for a file or a study, an association request is sent to the destination AE. Upon successful negotiation of a Presentation Context the object transfer is started and a C-STORE request is

 ROLAND CONSULT Electrophysiology and Imaging	DICOM Conformance Statement RETI-Port/Scan 21
No. 01 Rev. 01	

initiated for each object selected for storage. If the association cannot be opened, an error is displayed to the user. For each file sent in a send-job, the C-STORE response is received. If one or more C-STORE operations fail, a message is displayed to the user. For files that could not be sent because either the association could not be established or the C-STORE response status differs from the status SUCCESS or WARNING, or the transfer of the acquired documents fails for any other reason, the failed transfer attempt will be recorded in RETI-Port/Scan 21's database. A configurable number of retries will be executed in configurable time intervals until either the transfer succeeds or the maximum number of retries is reached.

The Storage process is initiated from the "Printout" dialog. "Send to PACS" is used instead of "Printout" is the operator handles an open worklist request. After printing all documents (PDF format) the operator calls "finish and send to PACS" from the DICOM window.

Error reporting: See 4.2.2

After all files of a send-job are sent, the association is released.

4.3.1.3.1.2 Proposed Presentation Contexts

The Send AE is capable of proposing the Presentation Contexts shown in the following table:

Table 4.3-6 PROPOSED PRESENTATION CONTEXTS FOR ACTIVITY SEND IMAGES

Presentation Context Table					
Abstract Syntax		Transfer Syntax		Role	Ext. Neg.
Name	UID	Name List	UID List		
Encapsulated PDF Storage	1.2.840.10008.5.1.4.1.1.104.1	Implicit VR Little Endian Explicit VR Little Endian	1.2.840.10008.1.2 1.2.840.10008.1.2.1	SCU	None

4.3.1.3.1.3 SOP Specific Conformance for Image Storage SOP Classes

All Image SOP Classes supported by the Send AE exhibit the same behavior, except where stated, and are described together in this section.

- 1 If none of the proposed Presentation Contexts is accepted by the SCP, then the Association is aborted by the Send AE using A-ABORT and the send job is marked as failed. The job failure is logged and reported to the user. The Send AE will retry failed transfers as described in 4.3.1.3.1.1.

Any of the C-Store response failures become reported to the user and the whole Store procedure may be repeated or the error shall be reported to the Roland Consult service desk.

The behavior of Send AE when encountering status codes in a C-STORE response is summarized in the Table below:

Table 4.3-7 STORAGE C-STORE RESPONSE STATUS HANDLING BEHAVIOR

Service Status	Further Meaning	Error Code	Behavior
Success	Success	0000	The SCP has successfully stored the SOP Instance. If all SOP Instances in a send job have status success then the job is marked as complete.

DICOM Conformance Statement

RETI-Port/Scan 21



No. 01
Rev. 01

Refused	Out of Resources	A700-A7FF	The Association is released using A-RELEASE-RQ and the send job is marked as failed. The status meaning is logged and the job failure is reported to the user. This is a transient failure.
Error	Data Set does not match SOP Class	A900-A9FF	The Association is released using A-RELEASE-RQ and the send job is marked as failed. The status meaning is logged and the job failure is reported to the user via the job control application.
Error	Cannot Understand	C000-CFFF	The Association is released using A-RELEASE-RQ and the send job is marked as failed. The status meaning is logged and the job failure is reported to the user.
Warning	Coercion of Data Elements	B000	Image transmission is considered successful but the status meaning is logged.
Warning	Data Set does not match SOP Class	B007	Image transmission is considered successful but the status meaning is logged.
Warning	Elements Discarded	B006	Image transmission is considered successful but the status meaning is logged.
*	*	Any other status code.	The Association is released using A-RELEASE-RQ and the send job is marked as failed. The status code is logged and the job failure is reported to the user.

The behavior of Storage AE during communication failure is summarized in the Table below:

Table 4.3-8 STORAGE COMMUNICATION FAILURE BEHAVIOR

Exception	Behavior
Timeout	The Association is released using A-RELEASE-RQ and the send job is marked as failed. The reason is logged and the job failure is reported to the user.
Association aborted by the SCP or network layers	The send job is marked as failed. The reason is logged and the job failure is reported to the user.

A failed send job is automatically restarted in periodic intervals until the maximum number of retries is reached.

The contents of all Image Storage SOP Instances created by RETI-Port/Scan 21 conform to the DICOM IOD definitions and are described in section 8.1.

4.3.2 Workflow Application Entity Specification

4.3.2.1 SOP Classes

The Workflow AE provides Standard Conformance to the following SOP Class:

 ROLAND CONSULT Electrophysiology and Imaging	DICOM Conformance Statement RETI-Port/Scan 21
No. 01 Rev. 01	

Table 4.3-9 SOP CLASSES FOR WORKFLOW AE

SOP Class Name	SOP Class UID	SCU	SCP
Modality Worklist Information Model – FIND	1.2.840.10008.5.1.4.31	Yes	No

4.3.2.2 Association Policies

4.3.2.2.1 General

The DICOM standard application context name for DICOM 3.0 is always proposed:

Table 4.3-10 DICOM APPLICATION CONTEXT FOR WORKFLOW AE

Application Context Name	1.2.840.10008.3.1.1.1
--------------------------	-----------------------

4.3.2.2.2 Number of Associations

The Workflow AE initiates one Association at a time for a Worklist request.

Table 4.3-11 NUMBER OF ASSOCIATIONS INITIATED FOR WORKFLOW AE

Maximum number of simultaneous Associations	1
---	---

The Workflow AE does not accept Association Requests by remote AEs.

4.3.2.2.3 Asynchronous Nature

The Workflow AE does not support asynchronous communication (multiple outstanding transactions over a single Association).

Table 4.3-12 ASYNCHRONOUS NATURE AS A SCU FOR WORKFLOW AE

Maximum number of outstanding asynchronous transactions	1
---	---

4.3.2.2.4 Implementation Identifying Information

The implementation information for this Application Entity is:

Table 4.3-13 DICOM IMPLEMENTATION CLASS AND VERSION FOR WORKFLOW AE

Implementation Class UID	1.2.826.0.1.3680043.2.891.113
Implementation Version Name	DICOM_CONNECT_25

4.3.2.1 Association Initiation Policy

The Workflow AE will initiate an Association in order to issue a C-FIND request according to the Modality Worklist Information Model

4.3.2.1.1 Activity – Worklist Update

4.3.2.1.1.1 Description and Sequencing of Activities

The request for a Worklist Update is initiated by user interaction. A user can select different search criteria, i.e. Scheduled Procedure Step Start Date, Modality and Scheduled Station AE Title (always set by default).

Parameters: Examination time range, patients PACS id (or * for all patients scheduled for this AET)

Upon initiation of the request, the Workflow AE will build a Request Identifier for the C-FIND request, initiate an Association to send the request with the search criteria and will wait for the Worklist responses. To protect the system from overflow, the Workflow AE can limit the number of processed worklist responses to a configurable maximum. By default, no limit is configured. Worklist items of the last request are temporarily stored and every worklist request updates the current list of open requests. So the requests can be transferred among the different RETI applications. During receiving the worklist response items are counted and the query processing is canceled by issuing a C-FIND-CANCEL if the configurable limit of items is reached. The results will be displayed in a list, which will be cleared with the next worklist update.





Figure 4-3 SEQUENCING OF ACTIVITY – WORKLIST UPDATE

A possible sequence of interactions between the Workflow AE and a Departmental Scheduler (e.g. a device such as a RIS or HIS which supports the Modality Worklist SOP Class as a SCP) is illustrated in the Figure above:

1. The Workflow AE opens an association with the Departmental Scheduler

 ROLAND CONSULT Electrophysiology and Imaging	<h1>DICOM Conformance Statement</h1> <p>RETI-Port/Scan 21</p>
No. 01 Rev. 01	

- The Workflow AE sends a C-FIND request to the Departmental Scheduler containing the Worklist Query attributes.
- The Departmental Scheduler returns a C-FIND response containing the requested attributes of the first matching Worklist Item.
- The Departmental Scheduler returns another C-FIND response containing the requested attributes of the second matching Worklist Item.
- The Departmental Scheduler returns another C-FIND response with status Success indicating that no further matching Worklist Items exist. This example assumes that only 2 Worklist items match the Worklist Query.
- The Workflow AE closes the association with the Departmental Scheduler.
- The user selects a Worklist Item from the Worklist and prepares to acquire new images.

4.3.2.1.1.2 Proposed Presentation Contexts

The Workflow AE will propose Presentation Contexts as shown in the following table:

Table 4.3-14 PROPOSED PRESENTATION CONTEXTS FOR ACTIVITY WORKLIST UPDATE

Presentation Context Table					
Abstract Syntax		Transfer Syntax		Role	Ext. Neg.
Name	UID	Name List	UID List		
Modality Worklist Information Model – FIND	1.2.840.10008.5.1.4.31	Implicit VR Little Endian Explicit VR Little Endian	1.2.840.10008.1.2 1.2.840.10008.1.2.1	SCU	None

4.3.2.1.1.3 SOP Specific Conformance for Modality Worklist

The Table below provides a description of the RETI-Port/Scan 21s Worklist Request Identifier and specifies the attributes that are copied into the images. Unexpected attributes returned in a C-FIND response are ignored.

Requested return attributes not supported by the SCP are set to have no value. Non-matching responses returned by the SCP due to unsupported optional matching keys are handled like matches of the request identifier. No attempt is made it filter out possible duplicate entries.

If an extended character set is used in the Request Identifier, Specific Character Set (0008,0005) will be included in the Identifier (see section 6 for supported values). Otherwise, Specific Character Set (0008,0005) will not be sent

The table below should be read as follows:

Attribute Name Attributes supported to build an RETI-Port/Scan 21 Worklist Request Identifier. RETI-Port/Scan 21 will supply this attribute as Return Key with zero length for Universal Matching, if it is not supplied as a matching key otherwise.

Tag DICOM tag for this attribute.

DICOM Conformance Statement

RETI-Port/Scan 21



No. 01
Rev. 01

- Match** Type of matching supported for this attribute by RETI-Port/Scan 21. Matching keys for Worklist Update. A "S" will indicate that RETI-Port/Scan 21 will supply an attribute value for Single Value Matching, a "R" will indicate Range Matching and a "*" will denote wild card matching.
- Query** An "x" will indicate that that RETI-Port/Scan 21 will supply this attribute as matching key automatically. An (x) will indicate that this matching key is provided only, if entered in the Query Patient Worklist dialog.
- Display** Displayed keys. An "x" indicates that this worklist attribute is displayed to the user during a patient registration dialog. For example, Patient Name will be displayed when registering the patient prior to an examination.
- IOD** An "x" indicates that this Worklist attribute is included into all Object Instances created during performance of the related Procedure Step.

Table 4.3-15: MODALITY WORKLIST - FIND SOP Class - C-FIND REQUEST IDENTIFIER

Attribute Name	Tag	Match	Query	Display	IOD	Note
Scheduled Procedure Step						
Scheduled Procedure Step Sequence	(0040,0100)					Only one item allowed
> Scheduled Station AE Title	(0040,0001)	S	x	x		
> Scheduled Procedure Step Start Date	(0040,0002)	S,R	(x)	x		
> Scheduled Procedure Step Start Time	(0040,0003)	S,R	(x)	x		
> Modality	(0008,0060)	S	(x)			
> Scheduled Physicians Name	(0040,0006)	S,*	(x)	x	x	Set as Performing Physician in IOD
> Scheduled Procedure Step Description	(0040,0007)			x		
> Scheduled Protocol Code Sequence	(0040,0008)					
>> Code Value	(0008,0100)					
>> Coding Scheme Designator	(0008,0102)					
>> Code Meaning	(0008,0104)			x		
> Scheduled Procedure Step ID	(0040,0009)					
> Scheduled Station Name	(0040,0010)	S	(x)	x		
> Scheduled Procedure Step Location	(0040,0011)	S	(x)	x		

No. 01
Rev. 01

> Pre-Medication	(0040,0012)			x		
> Requested Contrast Agent	(0032,1070)			x		
Requested Procedure						
Requested Procedure ID	(0040,1001)	S	(x)		x	Set as Study ID in IOD
Requested Procedure Description	(0032,1060)			x	x	Set as Study Description in IOD
Requested Procedure Code Sequence	(0032,1064)					Only one item allowed.
> Code Value	(0008,0100)					
> Coding Scheme Designator	(0008,0102)					
> Code Meaning	(0008,0104)			x		
Study Instance UID	(0020,000D)				x	
Patient Transport Arrangements	(0040,1004)			x		
Reason for the Requested Procedure	(0040,1002)					
Reason for Requested Procedure Code Sequence	(0040,100A)					
> Code Value	(0008,0100)					
> Coding Scheme Designator	(0008,0102)					
> Code Meaning	(0008,0104)					
Requested Procedure Priority	(0040,1003)					
Names Of Intended Recipients Of Results	(0040,1010)					
Requested Procedure Comments	(0040,1400)					
Imaging Service Request						
Accession Number	(0008,0050)	S	(x)	x	x	
Referring Physician's Name	(0008,0090)			x	x	
Requesting Physician	(0032,1032)			x		
Imaging Service Request Comments	(0040,2400)					
Visit Identification						
Admission ID	(0038,0010)				x	

DICOM Conformance Statement

RETI-Port/Scan 21



No. 01
Rev. 01

Institution Name	(0008,0080)				x	
Institution Address	(0008,0081)				x	
Institutional Department Name	(0008,1040)				x	
Visit Status						
Current Patient Location	(0038,0300)					
Patient Identification						
Patient's Name	(0010,0010)	S,*	(x)	x	x	
Patient ID	(0010,0020)	S	(x)	x	x	
Issuer of Patient ID	(0010,0021)	S	(x)	x	x	
Patient Demographic						
Patient's Birth Date	(0010,0030)			x	x	
Patient's Sex	(0010,0040)			x	x	
Patient's Size	(0010,1020)			x	x	
Patient's Weight	(0010,1030)			x	x	
Patient Medical						
Patient State	(0038,0500)			x		
Pregnancy Status	(0010,21C0)			x		
Medical Alerts	(0010,2000)			x		
Allergies	(0010,2110)			x		
Special Needs	(0038,0050)			x		
Patient Comments	(0010,4000)			x	x	
Ethnic Group	(0010,2160)			x	x	
Additional Patient History	(0010,21B0)			x		
Confidentiality Constraint on Patient Data Description	(0040,3001)			x		

The behavior of The Workflow AE, when encountering status codes in a Modality Worklist C-FIND response is summarized in the Table below. If any other SCP response status than "Success" or "Pending" is received by The Workflow AE, a message "query failed" will appear on the user interface.

 ROLAND CONSULT Electrophysiology and Imaging	<h1>DICOM Conformance Statement</h1> RETI-Port/Scan 21
No. 01 Rev. 01	

Table 4.3-16 MODALITY WORKLIST C-FIND RESPONSE STATUS HANDLING BEHAVIOR

Service Status	Further Meaning	Error Code	Behavior
Success	Matching is complete	0000	The SCP has completed the matches. Worklist items are available for display or further processing.
Refused	Out of Resources	A700	The Association is released using A-RELEASE-RQ and the worklist query is marked as failed. The status meaning is logged and reported to the user. Any additional error information in the Response will be logged.
Failed	Identifier does not match SOP Class	A900	The Association is released using A-RELEASE-RQ and the worklist query is marked as failed. The status meaning is logged and reported to the user. Any additional error information in the Response will be logged.
Failed	Unable to Process	C000 – CFFF	The Association is released using A-RELEASE-RQ and the worklist query is marked as failed. The status meaning is logged and reported to the user. Any additional error information in the Response will be logged.
Cancel	Matching terminated due to Cancel request	FE00	Worklist items are available for display or further processing. The status meaning is logged.
Pending	Matches are continuing	FF00	The worklist item contained in the Identifier is collected for later display or further processing.
*	*	Any other status code.	The Association is released using A-RELEASE-RQ and the worklist is marked as failed. The status meaning is logged and reported to the user. Any additional error information in the Response will be logged.

The behavior of The Workflow AE during communication failure is summarized in the Table below.

Table 4.3-17 MODALITY WORKLIST COMMUNICATION FAILURE BEHAVIOR

Exception	Behavior
Timeout	The Association is aborted using A-ABORT and the worklist query marked as failed. The reason is logged and reported to the user if an interactive query.
Association aborted by the SCP or network layers	The worklist query is marked as failed. The reason is logged and reported to the user if an interactive query.

Acquired images will always use the Study Instance UID specified for the Scheduled Procedure Step (if available). If an acquisition is unscheduled, a Study Instance UID will be generated locally.

4.3.2.1.1.4 Association Acceptance Policy

The Workflow AE does not accept Associations.

DICOM Conformance Statement

RETI-Port/Scan 21



No. 01
Rev. 01

4.3.3 Verification Application Entity

4.3.3.1 SOP Classes

The Verification Application Entity provides Standard Conformance to the following SOP Class:

SOP Class Name	SOP Class UID	SCU	SCP
Verification SOP Class	1.2.840.10008.1.1	Yes	Yes

4.3.3.2 Association Policies

The Verification AE can both accept and propose Association Requests for the Verification Service. The Verification AE will accept associations from any AET, no pre-configuration of known AEs is required. However, the Verification AE will reject associations if the called AE Title does not match the preconfigured AET of RETI-Port/Scan 21

4.3.3.2.1 General

The DICOM standard Application Context Name for DICOM 3.0 is always accepted and proposed

Table 4.3-18 DICOM APPLICATION CONTEXT FOR VERIFICATION AE

Application Context Name	1.2.840.10008.3.1.1.1
--------------------------	-----------------------

4.3.3.2.2 Number of Associations

The Verification AE initiates one Association at a time for each destination that DICOM connectivity shall be verified.

The association will be rejected if the maximum allowable number of associations is exceeded.

Table 4.3-19 NUMBER OF SIMULTANEOUS ASSOCIATIONS FOR VERIFICATION AE

Maximum number of simultaneous Associations accepted by the Verification AE	Unlimited (not configurable)
Maximum number of simultaneous Associations initiated by the Verification AE	1 (not configurable)

4.3.3.2.3 Asynchronous Nature

The Verification AE does not support asynchronous communication (multiple outstanding transactions over a single Association).

Table 4.3-20 ASYNCHRONOUS NATURE AS A SCU FOR VERIFICATION AE

Maximum number of outstanding asynchronous transactions	1
---	---

4.3.3.2.4 Implementation Identifying Information

 ROLAND CONSULT Electrophysiology and Imaging	DICOM Conformance Statement RETI-Port/Scan 21
No. 01 Rev. 01	

The implementation information for this Application Entity is:

Table 4.3-21 DICOM IMPLEMENTATION CLASS AND VERSION FOR VERIFICATION AE

Implementation Class UID	1.2.826.0.1.3680043.2.891.113
Implementation Version Name	DICOM_CONNECT_25

4.3.3.3 Association Initiation and Acceptance Policy

Due to the trivial functionality of the Verification SOP class, both, association initiation and acceptance policy are combined in this chapter. In Chapter 4.1.1 it is shown that “Verify connectivity with remote AE” is a symmetric activity which is executed by the local Verification AE and any remote AE in the same manner.

When the Verification AE accepts an association, it will respond to storage requests. The Verification AE will accept associations from any AET; no pre-configuration of known AEs is required. However, the Verification AE will reject associations if the Called AE Title does not match the preconfigured AET of the RETI-Port/Scan 21. The Verification AE may reject association requests for other reasons also. See 8.3 for details.

4.3.3.3.1 Activity: Verify connectivity with remote AE

4.3.3.3.1.1 Description and Sequencing of Activity

The Verification AE will initiate a new Association each time the user selects to verify connectivity for a specific AE.

Furthermore, it will respond to verification requests issued by a remote AE. Both scenarios are illustrated in Figure 4-4.




Figure 4-4 SEQUENCING OF ACTIVITY: VERIFY CONNECTIVITY WITH REMOTE AE

RETI-Port/Scan 21 verifies the connectivity to a Remote AE

1. The Verification AE opens an Association with the Remote AE

DICOM Conformance Statement

RETI-Port/Scan 21



No. 01
Rev. 01

2. The Verification AE sends a C-ECHO-Request to the remote AE. It awaits the C-ECHO-Response and propagates the result of the verification (successful/unsuccessful) to the user
3. The Verification AE closes the association

A Remote AE verifies the connectivity to RETI-Port/Scan 21

4. The Remote AE opens an Association with the Verification AE
5. The Remote AE sends a C-ECHO-Request to the Verification AE. The Verification AE will send a SUCCESSful response
6. The Remote AE closes the Association.

4.3.3.3.1.2 Proposed / Accepted Presentation Contexts

The Verification AE will propose and accept Presentation Contexts as shown in the following table:

Table 4.3-22 PROPOSED PRESENTATION CONTEXTS FOR VERIFICATION AE

Presentation Context Table					
Abstract Syntax		Transfer Syntax		Role	Ext. Neg.
Name	UID	Name List	UID List		
Verification	1.2.840.10008.1.1	Implicit VR Little Endian	1.2.840.10008.1.2	SCU / SCP	None
		Explicit VR Little Endian	1.2.840.10008.1.2.1		

4.3.3.3.1.3 SOP Specific Conformance for Verification SOP Class

The Verification AE provides standard conformance to the Verification SOP Class as an SCP and as an SCU. Technically, the Verification SCP is included in the system service process that is also hosting the Send AE. This means, that Verification requests can be handled as soon as the system service is started.

If the C-ECHO request was successfully received, a 0000 (Success) status code will be returned in the C-ECHO response. Otherwise, a C000 (Error - Cannot Understand) status code will be returned in the C-ECHO response.

In turn, the Verification AE treats a successfully received C-ECHO response with 0000 (Success) status code as a successful verification. Any other result will be treated as a verification failure.

4.4 NETWORK INTERFACES

4.4.1 Physical Network Interface

The DICOM Interface of the RETI-Port/Scan 21 provides DICOM TCP/IP Network Communication Support and uses the TCP/IP protocol stack from the operating system. It uses the Merge DICOM Toolkit subroutine library. All available Ethernet interfaces are supported.

4.4.2 Additional Protocols

Not applicable.

 ROLAND CONSULT Electrophysiology and Imaging	<h1>DICOM Conformance Statement</h1> <p>RETI-Port/Scan 21</p>
No. 01 Rev. 01	

5 CONFIGURATION

5.1 AE Title/Presentation Address Mapping

5.1.1 Local AE Titles

AE Title can be configured via the configuration interface. The default AE Title is RETI. All local AE's are using the same local AE Title.

Table 5.1-1 AE TITLE CONFIGURATION TABLE

Application Entity	Default AE Title	Default TCP/IP Port
Send	(must be configured)	Not Applicable
Storage Service		(must be configured)
Workflow		Not Applicable
Verification		(must be configured)

5.2 Parameters

A large number of parameters related to acquisition and general operation can be configured using the Service/Installation Tool or through editing configuration files manually. The Table below only shows those configuration parameters relevant to DICOM communication.

Table 5.2-1 CONFIGURATION PARAMETERS TABLE

Parameter	Configurable (Yes/No)	Default Value
General Parameters		
Max PDU Receive Size	Yes	64234 Bytes
Max PDU Send Size (larger PDUs will never be sent, even if the receiver supports a larger Max PDU Receive Size. If the receiver supports a smaller Max PDU Receive Size then the Max PDU Send Size will be reduced accordingly for the duration of the Association. Max PDU Receive Size information is exchanged during DICOM Association Negotiation in the Maximum Length Sub-Item of the A-ASSOCIATION-RQ and A-ASSOCIATE-AC)	Yes	64234 Bytes
Time-out waiting for association request or waiting for the peer to shut down an association. (ARTIM Timeout)	Yes	30 s
Time-out awaiting a reply to associate request	Yes	15 s
Time-out awaiting a reply to associate release	Yes	15 s
Time-out awaiting a network-write to be accepted.	Yes	15 s
Time-out awaiting a network-connect to be accepted.	Yes	15 s
Localization Parameters		
Specific Character Set used in messages	Yes	ISO_IR 6

DICOM Conformance Statement

RETI-Port/Scan 21



ROLAND CONSULT
Electrophysiology and Imaging

No. 01

Rev. 01

Parameter	Configurable (Yes/No)	Default Value
Retry of Failed Transfer Parameters		
Maximum number of DICOM objects transferred in a single association	Yes	100
Maximum number of transfer attempts until the transfer is considered as "failed"	Yes	5
Time Interval Between 2 Transfer Attempts for an SOP Instance	Yes	1 hour

No. 01**Rev. 01****6 SUPPORT OF CHARACTER SETS**

RETI-Port/Scan 21 support character sets

ISO_IR 6 (ISO 646 Default repertoire)ISO_IR 100 (ISO 8859-1 Latin Alphabet No. 1 supplementary set)

ISO_IR 101 (ISO 8859-2 Latin Alphabet No. 2 supplementary set)

ISO_IR 109 (ISO 8859-3 Latin Alphabet No. 3 supplementary set)

ISO_IR 110 (ISO 8859-4 Latin Alphabet No. 4 supplementary set)

ISO_IR 148 (ISO 8859-9 Latin Alphabet No. 5 supplementary set)

ISO_IR 144 (ISO 8859-5 Cyrillic Alphabet supplementary set)

ISO_IR 127 (ISO 8859-6 Arabic Alphabet supplementary set)

ISO_IR 126 (ISO 8859-7 Greek Alphabet supplementary set)

ISO_IR 138 (ISO 8859-8 Hebrew Alphabet supplementary set)

ISO_IR 13 (JIS X 0201: Katakana)

ISO_IR 166 (TIS 620-2533: Thai)

ISO_IR 192 (Unicode in UTF-8 for Chinese Alphabet)

GB18030

GBK

RETI-Port/Scan 21 does not support code extension techniques as described in ISO/IEC 2022:1994

DICOM Conformance Statement

RETI-Port/Scan 21



No. 01
Rev. 01

7 SECURITY

RETI-Port/Scan 21 does not support any specific security measures.

It is assumed that RETI-Port/Scan 21 is used within a secured environment. It is assumed that a secured environment includes at a minimum:

- a. Firewall or router protections to ensure that only approved external hosts have network access to RETI-Port/Scan 21.
- b. Firewall or router protections to ensure that RETI-Port/Scan 21 only has network access to approved external hosts and services.
- c. Any communication with external hosts and services outside the locally secured environment use appropriate secure network channels (e.g. such as a Virtual Private Network (VPN))

Other network security procedures such as automated intrusion detection may be appropriate in some environments. Additional security features may be established by the local security policy and are beyond the scope of this conformance statement.

8 ANNEXES

8.1 IOD CONTENTS

8.1.1 Created SOP Instances

The sections in this chapter specify the attributes of SOP Instances handled by RETI-Port/Scan 21. This refers to both, SOP Instances created by RETI-Port/Scan 21 and sent to remote AEs by the Send AE and SOP Instances received by the Storage Service AE to be processed by RETI-Port/Scan 21.

The following tables use a number of abbreviations. The abbreviations used in the “Presence of ...” column are:

VNAP	Value Not Always Present (attribute sent zero length if no value is present)
ANAP	Attribute Not Always Present
ALWAYS	Always Present
EMPTY	Attribute is sent without a value

The abbreviations used in the “Source” column:

MWL	the attribute value source Modality Worklist, see chapter 8.1.4
USER	the attribute value source is from User input
AUTO	the attribute value is generated automatically
CONFIG	the attribute value source is a configurable parameter

NOTE: All dates and times are encoded in the local configured calendar and time. Date, Time and Time zone are configured using the Service/Installation Tool.

8.1.1.1 Encapsulated PDF IOD

Table 8.1-1 IOD OF CREATED ENCAPSULATED PDF SOP INSTANCES

IE	Module	Reference	Presence of Module
Patient	Patient	Table 8.1-2	ALWAYS
	Clinical Trial Subject	-	NEVER
Study	General Study	Table 8.1-3	ALWAYS
	Patient Study	Table 8.1-4	ALWAYS
	Clinical Trial Study	-	NEVER
Series	Encapsulated Document Series	Table 8.1-9	ALWAYS
	Clinical Trial Series	-	NEVER

DICOM Conformance Statement

RETI-Port/Scan 21



No. 01
Rev. 01

IE	Module	Reference	Presence of Module
Equipment	General Equipment	Table 8.1-6	ALWAYS
	SC Equipment	Table 8.1-8	ALWAYS
Document	Encapsulated Document	Table 8.1-10	ALWAYS
	SOP Common	Table 8.1-7	ALWAYS

8.1.2 IOD Module Definitions

8.1.2.1 Common Modules

Table 8.1-2 PATIENT MODULE ATTRIBUTES OF CREATED SOP INSTANCES

Attribute Name	Tag	VR	Value	Presence of Value	Source
Patient's Name	(0010,0010)	PN	From Modality Worklist or user input. Values supplied via Modality Worklist will be entered as received. Maximum 64 characters.	VNAP	MWL/USER
Patient ID	(0010,0020)	LO	From Modality Worklist or user input. Maximum 64 characters.	VNAP	MWL/USER
Issuer of Patient ID	(0010,0021)	LO	From Modality Worklist or user input. Maximum 64 characters.	ANAP	MWL/USER
Patient's Birth Date	(0010,0030)	DA	From Modality Worklist or user input	VNAP	MWL/USER
Patient's Sex	(0010,0040)	CS	From Modality Worklist or user input	VNAP	MWL/USER
Patient Comments	(0010,4000)	LT	From Modality Worklist or user input	ANAP	MWL/USER

Table 8.1-3 GENERAL STUDY MODULE ATTRIBUTES OF CREATED SOP INSTANCES

Attribute Name	Tag	VR	Value	Presence of Value	Source
Study Instance UID	(0020,000D)	UI	From Modality Worklist or generated by device	ALWAYS	MWL/AUTO
Study Date	(0008,0020)	DA	Generated by device	ALWAYS	AUTO
Study Time	(0008,0030)	TM	Generated by device	ALWAYS	AUTO
Referring Physician's Name	(0008,0090)	PN	From Modality Worklist or user input	VNAP	MWL/USER
Study ID	(0020,0010)	SH	From Modality Worklist or user input	VNAP	USER

No. 01
Rev. 01

Accession Number	(0008,0050)	SH	From Modality Worklist or user input	VNAP	MWL/ USER
Study Description	(0008,1030)	LO	User Input	VNAP	USER

Table 8.1-4 PATIENT STUDY MODULE ATTRIBUTES OF CREATED SOP INSTANCES

Attribute Name	Tag	VR	Value	Presence of Value	Source
Admission ID	(0038,0010)	LO	From Modality Worklist or user input	VNAP	MWL/ USER

Table 8.1-5 GENERAL SERIES MODULE ATTRIBUTES OF CREATED SOP INSTANCES

Attribute Name	Tag	VR	Value	Presence of Value	Source
Modality	(0008,0060)	CS	DOC	ALWAYS	AUTO
Series Instance UID	(0020,000E)	UI	Generated by device	ALWAYS	AUTO
Series Number	(0020,0011)	IS	Generated by device	VNAP	AUTO
Performing Physician's Name	(0008,1050)	PN	Physician field in Study list. Maximum 64 characters.	ANAP	USER
Series Date	(0008,0021)	DA	Generated by device	VNAP	AUTO
Series Time	(0008,0031)	TM	Generated by device	VNAP	AUTO
Series Description	(0008,103E)	LO	User input	ANAP	USER

Table 8.1-6 GENERAL EQUIPMENT MODULE ATTRIBUTES OF CREATED SOP INSTANCES

Attribute Name	Tag	VR	Value	Presence of Value	Source
Manufacturer	(0008,0070)	LO	softgate	ALWAYS	AUTO
Institution Name	(0008,0080)	LO	From Configuration or Worklist	ANAP	CONFIG / MWL
Institution Address	(0008,0080)	ST	From Configuration or Worklist	ANAP	CONFIG / MWL
Station Name	(0008,1010)	SH	From Configuration	ANAP	CONFIG
Institutional Department Name	(0008,1040)	LO	From Configuration or Worklist	ANAP	CONFIG
Manufacturer's Model Name	(0008,1090)	LO	"RETI"	ANAP	AUTO
Software Version	(0018,1020)	LO	From Configuration, currently V10	ANAP	CONFIG

DICOM Conformance Statement

RETI-Port/Scan 21



No. 01
Rev. 01

Table 8.1-7 SOP COMMON MODULE ATTRIBUTES OF CREATED SOP INSTANCES

Attribute Name	Tag	VR	Value	Presence of Value	Source
Specific Character Set	(0008,0005)	CS	Configuration dependent, depending on the locale settings of the device. See Chapter 6.	ALWAYS	CONFIG
SOP Class UID	(0008,0016)	UI	Dependent on / according to the type of created object	ALWAYS	AUTO
SOP Instance UID	(0008,0018)	UI	Generated by device	ALWAYS	AUTO

8.1.2.1 Secondary Capture Image specific Modules

Table 8.1-8 SC EQUIPMENT MODULE ATTRIBUTES OF CREATED SOP INSTANCES

Attribute Name	Tag	VR	Value	Presence of Value	Source
Conversion Type	(0008,0064)	CS	WSD	ALWAYS	AUTO
Modality	(0008,0060)	CS	OT	ALWAYS	AUTO

8.1.2.2 Encapsulated PDF Specific Attributes

Table 8.1-9 ENCAPSULATED DOCUMENT SERIES MODULE ATTRIBUTES OF CREATED SOP INSTANCES

Attribute Name	Tag	VR	Value	Presence of Value	Source
Modality	(0008,0060)	CS	DOC	ALWAYS	AUTO
Series Instance UID	(0020,000E)	UI	Generated by device	ALWAYS	AUTO
Series Number	(0020,0011)	IS	Generated by device	VNAP	AUTO
Series Description	(0008,103E)	LO	User input	ANAP	USER
Protocol Name	(0018,1030)	LO	Generated by device	ANAP	USER


 ROLAND CONSULT Electrophysiology and Imaging	DICOM Conformance Statement RETI-Port/Scan 21
No. 01 Rev. 01	

Table 8.1-10 ENAPSULATED DOCUMENT MODULE ATTRIBUTES OF CREATED SOP INSTANCES

Attribute Name	Tag	VR	Value	Presence of Value	Source
Instance Number	(0020,0013)	IS	Generated by device	ALWAY	AUTO
Content Date	(0008,0023)	DA	Generated by device	ALWAY	AUTO
Content Time	(0008,0033)	TM	Generated by device	ALWAY	AUTO
Acquisition DateTime	(0008,002A)	DT	Generated by device	ALWAY	AUTO
Image Laterality	(0020,0062)	CS	User input	ANAP	USER
Burned in Annotation	(0028,0301)	CS	Yes (Header in printout PDF)	ALWAYS	USER
Document Title	(0042,0010)	ST	User input	VNAP	USER
Concept Name Code Sequence	(0040,E008)	SQ	User input	VNAP	USER

8.1.3 Used Fields in received IOD by application

Images received by the Storage Service AE are limited to the SOP Classes that can be handled by RETI-Port/Scan 21. Images received are referenced in the local database, which makes use of the conventional identification attributes to distinguish patients, studies, series and instances. The usage of attributes received via Modality Worklist is described in section 4.2.2.3.1.3.

8.1.4 Attribute mapping

The relationships between attributes received via Modality Worklist and stored in acquired images are summarized in Table 8.1-11. The format and conventions used in Table 8.1-11 are the same as the corresponding table in DICOM Part 4, Annex M.6 [DICOM].

Table 8.1-11 ATTRIBUTE MAPPING BETWEEN MODALITY WORKLIST AND CREATED INSTANCES

Modality Worklist	Instance IOD
Patient Name	Patient Name
Patient ID	Patient ID
Issuer of Patient ID	Issuer of Patient ID
Patient's Birth Date	Patient's Birth Date
Patient's Sex	Patient's Sex
Patient's Size	Patient's Size

DICOM Conformance Statement

RETI-Port/Scan 21



No. 01
Rev. 01

Patient's Weight	Patient's Weight
Ethnic Group	Ethnic Group
Patient Comments	Patient Comments
----	----
----	Series Instance UID
----	Series Description
----	Protocol Name
----	Operator's Name
----	----
----	SOP Class UID
----	SOP Instance UID
Scheduled Performing Physician's Name	Performing Physician's Name
Referring Physician's Name	Referring Physician's Name
Admission ID	Admission ID
----	----
Study Instance UID	Study Instance UID
Accession Number	Accession Number
Scheduled Protocol Code Sequence	----
Requested Procedure Description	Study Description
Requested Procedure ID	Study ID
Scheduled Procedure Step ID	----
Modality	Modality ²
Scheduled Station AE Title	----
Scheduled Station Name	----
Scheduled Procedure Step Location	----
Scheduled Procedure Step Description	----
Scheduled Protocol Code Sequence	----
Requested Procedure Code Sequence	----
Institution Name	Institution Name
Institution Address	Institution Address
Institutional Department Name	Institutional Department Name

8.1.5 Coerced/Modified Fields

No coercion / modification of fields received by other DICOM AEs is performed by the AEs which are part of RETI-Port/Scan 21.

8.2 PRIVATE TRANSFER SYNTAXES

No Private Transfer Syntaxes are supported.

DICOM Conformance Statement

RETI-Port/Scan 21



No. 01
Rev. 01

8.3 Association Rejection Reasons

The Verification AE may reject association requests as shown in Table 8.3-1. The Result, Source and Reason / Diag columns represent the values returned in the appropriate fields of an ASSOCIATE-RJ PDU (see Section 9.3.4 "A-ASSOCIATE-RJ PDU Structure" in PS3.8). The contents of the Source column are abbreviated to save space and the meanings of the abbreviations are:

- a. 1 - DICOM UL service-user
- b. 2 - DICOM UL service-provider (ASCE related function)
- c. 3 - DICOM UL service-provider (Presentation related function)

Table 8.3-1 ASSOCIATION REJECTION REASONS

Result	Source	Reason / Diag	Explanation
2 – rejected transient	C	2 – local limit exceeded	The (configurable) maximum number of simultaneous association has been reached. An association request with the same parameters may succeed at a later time
2 – rejected permanent	A	2 – application context name not supported	The association request contained an unsupported Application Context Name. An association request with the same parameters will not succeed at a later time

9 Revision History

Rev	gültig ab	gültig bis	Änderung
01	15.03.2021	-	Release